Instructions for Use

Single Use Only | Disposable

Product Description

Hu-Friedy Bonded Double Wraps are square or rectangular nonwoven sheets produced using a three-layer SMS (spunbond-meltblown-spunbond) process. It consists of two sheets of SMS wrap, ultrasonically bonded together along two edges for convenient simultaneous wrapping of one or a collection of medical devices that will be sterilized following standard healthcare practices. They are composed of polypropylene with blue or green pigment and an anti-static treatment. The Bonded Double Wrap allows a sterilized package of medical devices to be opened aseptically, and supplied in a variety of sheet sizes as detailed in Table 1.

Table 1: Available Sizes of Hu-Friedy Bonded Double Wrap

Part Code	Dimensions	200 Grade
IMS-2215	15 x 15" (380 x 380mm)	X
IMS-2220	20 x 20" (508 x 508mm)	X
IMS-2224	24 x 24" (610 x 610mm)	X

Indications For Use

Hu-Friedy Bonded Double Wraps are intended to be used to enclose another medical device that is to be sterilized by a healthcare provider via the following:

- Pre-vacuum Steam 270°F/132°C for 4 minutes
 - Validated for dry time of 20 minutes for the 200 grades.
 - Validated for pre-vacuum steam sterilization of two stainless steel lumens 3 mm in diameter or larger and 400 mm in length or less.
- Gravity Steam 250°F/121°C for 30 minutes
 - Validated for dry times of 20 minutes.
 - Validated for gravity steam sterilization of two stainless steel lumens 3 mm in diameter or larger and 400 mm in length or less.
- 100% Ethylene Oxide (EO) with a concentration of 725-735 mg/L @ 131°F/55°C and 40% - 80% relative humidity for 60 minutes
 - Validated for aeration time of 8 hours at 55°C.
 - Validated for ethylene oxide sterilization of two stainless steel lumens 3 mm in diameter or larger and 400 mm in length or less.
- Advanced Sterilization Products (ASP) STERRAD® 100S
 - Validated with ASP STERRAD® 100S sterilization of lumens 2 mm inside diameter or larger and a length of 250mm or less.
- Advanced Sterilization Products (ASP) STERRAD® 100NX (Standard, Express and Flex cycles)
 - Validated with ASP STERRAD® 100NX sterilization cycles detailed in Table 2.
- STERIS Amsco® V-PRO 1, STERIS Amsco® V-PRO 1 Plus, STERIS Amsco® V-PRO maX Low Temperature Sterilization Systems (Lumen, Non Lumen and Flexible Cycles)
 - Validated for use with STERIS V-PRO® cycles detailed in Table 3.
 - Validated to be effectively aerated during the pre-programmed STERIS Amsco V-PRO® Sterilization Cycles.

The wrap is intended to allow sterilization of the enclosed medical device(s) and maintain sterility of the enclosed device(s) until used.

Table 2: Validated Advanced Sterilization Products (ASP) STERRAD® 100NX Cycles

Advanced Sterilization Products (ASP) STERRAD® System and Cycle	Intended Load
100NX Standard Cycle	Reusable metal and non-metal medical devices, including up to 10 lumens of the following dimensions per chamber load: <ul style="list-style-type: none"> • An inside diameter of 0.7mm or larger and a length of 500mm or shorter of single-channel stainless steel lumens.
100NX Flex Cycle	One or two single-channel flexible endoscopes with or without a silicone mat. Flexible endoscope may contain: <ul style="list-style-type: none"> • A single channel Teflon/PE/PTFE lumen with an inside diameter of 1mm or larger and a length of 850mm or shorter. (A maximum of two flexible endoscopes, one per tray per sterilization cycle)
100NX Express Cycle	Non-lumened reusable metal and non-metal medical devices requiring surface sterilization, or sterilization of mated stainless steel and titanium surfaces, and rigid or semi-rigid endoscopes without lumens.

Table 3: Validated STERIS Amsco® V-PRO Cycles

STERIS Amsco® V-PRO Cycle	Intended Load
Lumen Cycle	Reusable metal and non-metal medical devices, including up to 20 stainless steel lumens with dimensions of 3.0mm diameter or larger and a length of 400mm or shorter.
Non Lumen Cycle	Non Lumened reusable metal and non-metal medical devices.
Flexible Cycle	Single lumen surgical flexible endoscopes and bronchoscopes in the following load configuration: <ol style="list-style-type: none"> 1.) One tray containing a flexible endoscope with a light cord (if not integral to endoscope) and mat. The flexible endoscope having an inside diameter of 1mm or larger and a length of 8.50mm or shorter. 2.) An additional tray containing non-lumened medical devices.

Warnings:

- Do not apply wrap to the following sterilization methods: 1) Dry heat; 2) Radiation
- Inspect wrap for damage or extraneous matter prior to use. If a defect is detected, do not use.
- Do not use wrapped contents if wrap is torn or wet.

Manufactured for:
Hu-Friedy Mfg. Co., LLC
3232 North Rockwell St.
Chicago, IL 60618
USA

US Customer Care
1-800-HUFRIEDY
(1-800-483-7433)
Website: Hu-Friedy.com

Precautions:

- Do not open carton or package with a sharp knife, knives can easily damage the wrap.
- Prior to use, assure that all medical devices intended to be sterilized while wrapping within Hu-Friedy Bonded Double Wrap are compatible with and sterilizable by the sterilization modality and cycle listed in the Indications for Use within these instructions. Consult the sterilization instructions for all devices intended for sterilization. Some medical devices, regardless of the sterilization method and sterilization wrap/container applied, may require special considerations in packing configurations to ensure sterilization. Some of these considerations are referenced in ANSI/AAMI ST79 Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities.
- Do not use in the presence of flammable anesthesia.
- If sterilization is performed by an outside contract facility, Hu-Friedy recommends that the wrapped devices be protected from contamination by an additional covering.

Directions for Use

Hu-Friedy Bonded Double Wrap should be used in accordance with the preparation, wrapping and sterilization chamber loading recommendations of the following standards:

- ANSI/AAMI ST79: Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities
- CDC Guidelines for Infection Control in Dental Health-Care Settings – 2003
- ANSI/AAMI ST41: Ethylene Oxide Sterilization in Health Care Facilities
- AORN Standards, Recommended Practices, and Guidelines

Storage Prior to Use

- Storage location should meet the following conditions: 1) Clean; 2) Free of dust; 3) Away from fluorescent or ultraviolet light.

Prior to Use

- Recommended temperature and humidity conditions: 1) Minimum of two hours; 2) Temperature (68°F to 73°F) / (20°C to 23°C); 3) Relative humidity ranging from 30% to 60%
- Examine wrap and discard if damage or extraneous matter is detected.
- Thoroughly clean and dry items to be wrapped.

Wrapping with Hu-Friedy Bonded Double Wrap:

- Place item(s) on the double-sheet and wrap using typical aseptic wrapping technique per ANSI/AAMI ST79 and CDC Guidelines. Recommendations for wrap contents are provided in Table 4.
- Ensure that the first fold is pulled far enough to cover all package contents. Caution: Covering all package contents with the first fold is required for sterility maintenance, and failure to follow this for correct wrapping technique could compromise sterility.
- Securing the Wrapped Package:
 - Secure with an appropriate closure (tape or alternate closure suitable for the sterilization method to be applied) and label.
 - Closure must: 1) Allow the sterilant to penetrate the wrapped package; 2) Avoid constriction of the package; 3) Maintain package integrity.

Table 4: Wrap Model Recommendations²

Sterilization Wrap Models	Intended Load	Pre-Vacuum and Gravity* Steam and Ethylene Oxide	Advanced Sterilization Products (ASP) Sterrad® 100S	Advanced Sterilization Products (ASP) Sterrad® 100NX	Steris Amsco® V-PRO Lumen and Non Lumen Cycle	Steris Amsco® V-PRO Flex Cycle
200 Grade	Light Weight Package	6 lbs	10.7 lbs	10.7 lbs	10.7 lbs	6.5 lbs

Sterilization

- Hu-Friedy Bonded Double Wrap intended for use with the common healthcare sterilization parameters described in the Indications for Use. Users should consult the sterilizer manufacturer for appropriate sterilizer loading configurations.
- If a sterilizer malfunctions or the full cycle is not completed for any reason, packages should be re-wrapped with an unprocessed wrap and the sterilization to be repeated.
- Note: Many factors can affect drying time other than sterilization wrap, such as the pack configuration that is used, sterilizer loading, cycle variations, the performance of the sterilizer, temperature and steam distribution, altitude, and ambient environmental conditions. Sterilizers vary widely in design and performance characteristics. User should consult the sterilizer manufacturer's operator manual for specific drying times.

Post-Sterilization Cooling/Unloading

- Leave wrapped packages in the sterilizer until cool to avoid compromising package sterility.
- Visually inspect wrapped items as they are removed from the cart. Items that are damaged or wet should not be used.

Sterility Maintenance

- Real-time testing supports maintenance of package sterility for 180 days provided package integrity is maintained.
- Store wrapped packages as recommended in the ANSI/AAMI.
- Refer to and follow "Storage Prior to Use" and "Prior to Use" recommendations.
- Caution: Do not stack trays. Stacking trays can result in damage to the wrap caused by undue pressure from the weight.

Opening

- Inspect package for damage, wetness, or any sign of potential contamination prior to opening and again after opening but before use of contents. Caution: Do not use contents if these conditions are present, as sterility could be compromised. Reprocess the contents using an unprocessed wrap if any of these conditions are noted.
- Open packages aseptically in accordance with the healthcare facility's policy.

Disposal

- Do not re-use Hu-Friedy Bonded Double Wrap. Hu Friedy does not endorse the re-use (re-sterilization) of its sterilization wrap and does not warrant performance if product is re-used.
- Recycle, landfill or incinerate used wrap based upon state and local regulations. Note the wrap is composed of polypropylene plastic which has a plastics recycle code of "5".

Footnotes

- ¹ It is recommended the wrap layer containing the printed product information be visible as the outer layer after wrapping.
² Individual results may differ due to factors such as variations in wrapping techniques, handling practices and folding methods. Results may also differ due to the use of irregularly shaped contents that may result in added stress on the wrap. Each healthcare facility should determine for itself which wrap model is most appropriate for each intended use. ³ It is recommended to not exceed the maximum wrapped package content weights indicated for each wrap model. In addition, it is recommended to not exceed the number, weight and size of individual content types that were validated.

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